Mesa Italia S.r.L. via dell'Artigianato, 35/37, 25039 Travagliato (BS)
ITALY Tel. +39 0306863251 Fax +39 0306863252
info@mesaitalia.it sales@mesaitalia.it www.mesaitalia.it

MAGNUM SOLARE CE 0123

LOT00000

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MAGNUM SOLARE CE 0123

LOT 00000

Cobalt based dental alloy for ceramic, type 4*

*According to ISO 22674:2016

The product is intended for dental laboratories for the creation of products that must be rigid and resistant, single crowns, bridges and superstructures on implants.

The alloy MAGNUM SOLARE is distinguished by an outstanding fluidity, which grants to fill even the thinnest details of the frame, down to three tenth of millimeter.

Their molecular structure allows you to obtain smooth, compact surfaces with little oxide formation even after repeated oven treatment and their thermal expansion coefficient is ideal for every ceramic of the latest generation.

The alloy MAGNUM SOLARE is produced in conformity with standards ISO 9693:2020 and ISO 22674:2016, it is highly corrosion resistant, it does not contain Nickel and it is completely free of toxic elements like beryllium, indium and gallium. MAGNUM SOLARE is also supplied in discs for CAD/CAM manufacturing.

MAGNUM SOLARE is a product intended to be used by qualified technician.

Instruction for use

Modeling

Modeling with a minimum thickness of 0.3mm on simple structures and 0.5mm in complex cases or with patients with bruxism. Avoid elbows.

Casting Sprues

Direct pinning: Pin: Ø 3mm, Nourice: Ø >= 6mm

Indirect pinning with the stabilizer bar: Pin: Ø 3mm, Bar >= 5mm, Casting sprue >= 6mm

Casting

Cast MAGNUM SOLARE alloys in pre-heated and cleaned ceramic crucibles used for this alloy only. Each alloy should have its crucible; do not overheat the alloy and do not use the flux.

Reach a temperature of 850-950°C; maintenance time in pre-heated oven depends on investment characteristics and on the dimension of the cylinder.

Induction casting: start the casting when the cylinders will melt together and just before the superficial crust opens.

Open flame melting: place ingots in the pre-heated ceramic crucible and heat them evenly with circular movements. When ingots have melted, start centrifugal unit. Use multi-flame welding torches only. Do not use any flux. Indicative values for flame regulation: Acetylene 0,4 bar / Oxygen 2 bar; Propane 0,2 bar / Oxygen 2 bar.

Manufacturing

Let the cylinder cool down at room temperature, remove the investment and sandblast with aluminium oxide 110 to 150 µm at a pressure of 3-4 bar.

It is recommended an homogenization and degassing treatment of the casted piece at 1000°C for 10 minutes in atmosphere.

Continue manufacturing process using tungsten carbide burs.

The worked framework has to be sandblasted with disposable equipment in aluminium oxide of 110 to 150 µm at a max. pressure of 3-4 bar and then steam clean. After cleaning, the framework should not be touched anymore and should be held with hemostatic forceps only.

Oxidation

Fire for 5 min. under vacuum at 950-980°C. After firing, the oxide layer has to be carefully sandblasted using disposable oxide-aluminium equipment 110 to 150 µm at a pressure of 2.5-4 bar. Steam and check that the surface has an homogeneous grey surface.

Re-use of sprues

The best results are obtained with pure MAGNUM SOLARE. It is nonetheless possible to re-use sprues and cones only once, given that new metal coming from the same lot is added in equal quantity.

Ceramization

It's possible to use Bonding. Apply the ceramic, following its manufacturer's instructions. Slow cooling is recommended.

Soldering

We suggest using our MAGNUM SALDATURA Co to weld.

Waste disposal

The processing scrap must be disposed of as special waste in accordance with the EC directives 2008/98/CEE on waste, and 94/62/CEE on packaging and packaging waste and in compliance with national legislation in force on the subject.

Safety instructions

- Cobalt-Chrome-based alloys can seldom cause dermatitis on sensitive subjects. A Patch-Test is thus advisable.
- Metal dusts and smoke are dangerous for health. Use exhaust fans while casting and polishing.
- Before prosthesis application verify if other metal implants are in patient's oral cavity. Coexistence of different metals can cause a "pile" effect.
- We identify every batch with a number. We recommend to write it down in patient's file to allow its complete traceability.
- This product does not need any special preservation precautions. MESA keeps batch traceability for 15 years.
- The methods of preserving the products made with MAGNUM SOLARE are under the responsibility of professional operators (dental technicians). This product must be stored, transported and used at a temperature between -30°C and +50°C to ensure the integrity of the packaging and the instructions contained therein.
- MESA recommends that the patient be made aware of the possibility for dental alloys to affect MRI results.
- The product is sold non-sterile.

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